

TROM Adjuster™

POST-OPERATIVE BRACE

DONJOY®



TROM Adjuster™

POST-OPERATIVE BRACE

Indications

- Meniscal repairs
- Meniscal transplants
- Regenerative Chondroplasty
- Marrow stimulation/microfracture
- Osteochondral Grafting (allograft and autograft)
- Autologous chondrocyte implantation (ACI)
- Osteotomy (DFO/HTO)

Features & Benefits

- Off-loading adjustable hinge technology
- Tele-Fit: Sliding mechanism allowing independent strap movement just proximal and distal of the knee
- Telescoping adjustability between 17 3/4" and 24"
- Accommodates cold therapy wraps/dressings
- Sili grip strap pads assists in increasing brace suspension
- Warranty: 6 months

Suggested HCPCS code - L1832



Adjust. Off-load. Heal quickly. Heal safely.

Create the optimal knee off-loading with just a few turns. TROM Adjuster, another first-to-market product from DonJoy, offers controlled and consistent off-loading to create a better healing environment. Using our successful TROM platform, TROM Adjuster represents a comprehensive post-operative solution for meniscal and chondral defect procedures.



Off-loading adjustable hinge technology



+ adjustment
(Up to 9mm or ~ 3° of hinge adjustability)



- adjustment



Sili grip strap pads

ORDERING INFORMATION

Part Number	Description
11-5115-9	TROM Adjuster,Medial RT/Lateral LT
11-5116-9	TROM Adjuster,Medial LT/Lateral RT
11-5117-9	TROM Adjuster XL,Medial RT/Lateral LT
11-5118-9	TROM Adjuster XL,Medial LT/Lateral RT
11-1415	Sili grip strap pad kit (26" length)
25-2025	Replacement strap kit
25-2026	Replacement strap kit, XL
11-2024	Foam Liner Pad kit



DJO, LLC • 1430 Decision Street
Vista, CA 92081-8553 U.S.A.

Phone: 760.727.1280 • Toll Free: 800.336.6569
www.donjoy.com

