

ARMOR

ARMOR is the strongest, off-the-shelf functional knee brace on the market and is designed for active people seeking ligament protection, excellent fit, and superb comfort.

Ideal Applications

- Active adolescents and adults
- Contact, extreme, and water sports
- Moderate to severe ACL or PCL instabilities
- ACL/PCL combined instabilities (CI)
- ACL or PCL reconstructions
- MCL/LCL instabilities
- Hyperextension
- Prophylactic use

Features and Benefits

- Four-Points-of-Leverage design pioneered by dj Orthopedics
- FourcePoint™ Hinge Technology provides an adjustable, gradual “soft” stop to extension
- Standard polycentric hinge available on all versions
- Flexion and extension control
- Medial hinge deflector (Ski version only) for bilateral brace use
- Moldable frame for custom fit
- Medial swooping thigh cuff for inner thigh clearance and bilateral brace use
- Strong and lightweight 6061-T6 aircraft aluminum frame for medial / lateral control
- Durable, silver-veined powder coat finish
- Standard calf length for maximum leverage and short calf length for boot clearance
- Swiveling strap tabs to accommodate changing leg and muscle positions
- Anti-migration Supra Condyle Pad (Action version only)
- Air Condyle Pad (Ski version only) promotes brace suspension
- Impact Guard accessory protects the shin and knee against blunt blows
- Calf Pinch Guard accessory provides a smooth transition between boot and calf cuff
- 1 year warranty on frame and hinge
- 6 month warranty on soft goods



Ordering Information

PART NUMBER	LIGAMENT	HINGE	CALF LENGTH	LEG	SIZES AVAILABLE: MEASURE 6" ABOVE MID-PATELLA		
Armor Action					x = 1 = XS	13-15.5"	(33-39 cm)
NEW 11-1440-x	ACL	FourcePoint™	Short	Right	x = 2 = S	15.5-18.5"	(39-47 cm)
NEW 11-1441-x	ACL	FourcePoint™	Short	Left	x = 3 = M	18.5-21"	(47-53 cm)
NEW 11-1442-x	ACL	FourcePoint™	Standard	Right	x = 4 = L	21-23.5"	(53-60 cm)
NEW 11-1443-x	ACL	FourcePoint™	Standard	Left	x = 5 = XL	23.5-26.5"	(60-67 cm)
11-1004-x-26000	ACL	Standard	Short	Right	x = 6 = XXL	26.5-29.5"	(67-75 cm)
11-1005-x-26000	ACL	Standard	Short	Left	x = 7 = XXXL	29.5-32"	(75-81 cm)
NEW 11-1444-x	ACL	Standard	Standard	Right			
NEW 11-1445-x	ACL	Standard	Standard	Left			
11-1025-x-26000	PCL	Standard	Short	Right			
11-1026-x-26000	PCL	Standard	Short	Left			
NEW 11-1449-x	PCL	Standard	Standard	Right			
NEW 11-1450-x	PCL	Standard	Standard	Left			
11-1029-x-26000	CI	Standard	Short	Right			
11-1030-x-26000	CI	Standard	Short	Left			
NEW 11-1447-x	CI	Standard	Standard	Right			
NEW 11-1448-x	CI	Standard	Standard	Left			
Armor Ski							
NEW 11-1451-x	ACL	FourcePoint™	Short	Right			
NEW 11-1452-x	ACL	FourcePoint™	Short	Left			
11-9096-x-26000	ACL	Standard	Short	Right			
11-9097-x-26000	ACL	Standard	Short	Left			
11-1023-x-26000	PCL	Standard	Short	Right			
11-1024-x-26000	PCL	Standard	Short	Left			
11-1027-x-26000	CI	Standard	Short	Right			
11-1028-x-26000	CI	Standard	Short	Left			
					ACCESSORIES	DESCRIPTION	SIZE
					11-0075-x-00000	Black Lycra® Undergarment	S-XXL
					11-0122-x-00000	Black Neoprene Undergarment	XS-XXL
					11-1006-x-06000	Impact Guard	S-XXL
					11-1095-9-06000	Calf Pinch Guard	Universal
					11-0956-0-06000	Supra Condylar System	Universal
					11-0398-9-00000	Full Extension Stop	Universal
					12-0422-0-00000	Flexion/Extension Stop Kit	Universal
					25-2428	Mesh Brace Bag	Universal