

HAT TRICK™

HEADGUARD

DONJOY®



HAT TRICK™

HEADGUARD

Patients

- Children, Adolescents, & Adults

Indications

- Head-to-head impact protection while playing soccer

Features and Benefits

- Meets every performance requirement in ASTM F2439, "Standard Specification for Head gear Used in Soccer"
- Permitted by FIFA, the U.S. Soccer Federation, and National Fed. of State High School Associations
- Impact-absorbing foam layer covers and protects the forehead and temples without altering ball rebound speeds
- CoolMax® liner fabric accelerates moisture evaporation to keep the head cool
- Lycra outer fabric improves conformity to natural head shape
- 4 adjustable, latex-free Velcro straps maximize comfort and provide a uniform fit in the occipital area
- Ultra-lightweight (1.6 oz.), low profile, ventilated design allows for unrestricted, natural head and neck movement
- 6 month warranty



DONJOY

A State of Mind.

The DonJoy HAT TRICK headguard is clinically proven to reduce the risk of concussion resulting from head-to-head impact and is ideally suited for soccer players of all ages and ability levels.

HAT TRICK




ORDERING INFORMATION

Part Number	Size
11-1554	Small / Medium
11-1555	Large

Available in Black only.

SIZING



Size	Circumference Just Above Ears / Eyebrows	Cap Size
Small / Medium	19 1/4" – 23" (49 – 56 cm)	6 1/8 – 7
Large	22 3/4" – 24 1/2" (55 – 62 cm)	6 7/8 – 7 3/4



1430 Decision Street, Vista, CA 92081-8553
 t800.336.6569 f800.936.6569
 www.donjoy.com © 2006 DJO, LLC



It is within the sole discretion of the customer to determine the appropriate billing code for a product, as well as, whether the use of a product complies with medical necessity and other documentation requirements of the payor. DJO, LLC accepts no responsibility whatsoever in this regard. Withnall, C. et al. "Effectiveness of Headgear in Football." *British Journal of Sports Medicine*. 39: 40-48, 2005. DonJoy is a registered trademark and HAT TRICK is a trademark of DJO, LLC. CoolMax is a registered trademark of INVISTA.