

CCMI Mark II

Ideal Application

The Custom Contour Measuring Instrument is used when a Custom Defiance or a Custom Monarch Osteoarthritis brace is recommended. The device takes into account fifteen separate anatomical reference points from the patient's leg indicating knee joint angle, thigh and calf circumference, thigh and calf contours as well as knee width and varus/valgus

Features and Benefits

- Accurate and faster fabrication process
- Standard and short measuring options
- Precise soft-tissue measurement
- Easy and reliable to use
- Time efficient
- No mess to clean up, no use of casting products
- Own carrying case for easy storage
- No mailing of cast, just fax form in

Ordering Information

PART NUMBER

11-0192-9-06000

