

Ideal Application

This brace is ideal for patients recovering from injury or surgery, providing full and limited range-of-motion settings. With its off-set hinges designed to closely follow true knee kinematics, early mobilization can be facilitated through the convenient dial-in adjustment of flexion or extension settings.

Features and Benefits

- Patented swing-cap dual axis hinge
- Full and limited range-of-motion settings at 0°, 15°, 30°, 60° and 90°
- Universal sizing
- Integrated strap/hinge design for circumferential compression
- Breakdown capability on long version

Ordering Information

PART NUMBER	SIZE	LENGTH
11-0523-2-06060	Short/Univ	18" (45.75 cm)
11-0523-4-06060	Long/Univ	23" (58.50 cm)

Short = up to 28" inseam. Long = greater than 28" inseam.
Fits to 30" circumference

