

Contoured Wrist/Forearm

Ideal Application

The Contoured Wrist/Forearm Support provides comfort for post-surgical immobilization, post-fracture casting, and wrist sprains and strains.

Features and Benefits

- Vinyl construction
- Malleable dorsal and volar stays
- Pressure-sensitive closures for easy application

Ordering Information

PART NUMBER	
11-0218-x	Right
11-0219-x	Left
SIZES AVAILABLE: XS-XL PALM WIDTH SIZING	
x=1=Xsmall=2-2 3/4" (5-7 cm)	x=2=Small=2 3/4-3" (7-7.5 cm)
x=3=Medium=3-3 3/4" (7.5-9.5 cm)	x=4=Large=3 3/4-4 1/4" (9.5-10.75 cm)
x=5=XLarge=4 1/4-4 3/4" (10.75-11.25 cm)	
Color Available: grey	

