

## Telescoping Cool ELS®

### Ideal Application

For increased patient comfort and compliance in hot and humid climates, the ELS features a slide lock mechanism providing the ability to lock at full extension or have free range of motion. The brace is very easy to use with a simple yet effective design. The brace is ideal for post-op use when full extension weight-bearing protocols are used. The Cool ELS gives improved patient comfort with abbreviated foam.

### Features and Benefits

- Lock at full extension
- Easily size and custom fit the brace to the patient
- Brace can extend to ankle to help prevent migration
- 50 degree flexion stop kit available for meniscal repair protocols

### Ordering Information

*PART NUMBER*

11-0967-9-13066	Enhanced ELS, Cool
-----------------	--------------------

*ACCESSORIES*

11-0080-9-00000	50° flexion stop, ELS
-----------------	-----------------------

11-0024-x	Ankle stirrup, right
-----------	----------------------

11-0026-x	Ankle stirrup, left
-----------	---------------------

