

Defiance® Ideal Application

The Custom Defiance is DonJoy's hallmark knee brace. Combining the ultimate in technology, materials and design, this lightweight brace provides durable support for moderate to severe ACL, PCL, MCL and LCL instabilities. Plus, the low-profile design fits comfortably underneath clothing and sports uniforms, including baseball and football. As the strongest of DonJoy's functional braces, the custom, carbon composite Defiance is recommended for high collision sports, while its featherweight materials also make it the brace of choice for jumping sports, when every ounce counts. Available in three lengths; standard, short and standard thigh with short calf.

Features and Benefits

- ACL/PCL versions with patented 4-point dynamic leverage system
- Weighs approximately 18 ounces
- Low-profile design
- Custom made to individual patient measurements
- 72-hour in-house turnaround
- CI version with anterior thigh strap
- Short calf version available for ski and motocross boot clearance
- Short and standard lengths
- Defiance Extra version includes reinforced carbon frame and slide shield
- Suede, chamois and chamois pneumatic liners and strap pads available
- 16 individual color options, 6 graphic patterns, 18 graduation options
- Flexion stops: 45°, 60°, 75° and 90°
- Extension stops: 0°, 10°, 20°, 30° and 40° (10° installed)
- Easy to use custom measuring system (CCMI)
- Lifetime limited warranty on frame and hinge

Ordering Information

EACH DEFIANCE IS CUSTOM MADE.

| SIZE | THIGH | CALF |
|----------|------------------|------------------|
| Minimum: | 12 1/2" (31.8cm) | 10 1/4" (26cm) |
| Maximum: | 30 5/8" (77.8cm) | 20 3/4" (52.7cm) |

Measurements taken with CCMI Mark II or casted in weight bearing

Colors Available: red, yellow, orange, blue, metallic red, metallic dark blue, light blue, green, white, black, metallic bright blue, metallic green, metallic gold, metallic black, metallic purple, purple, red-white-blue, black with any flat color, white with any flat color, snake skin, leopard blue, tiger, tiger red, splash, and butterfly (see chart page 119)

ACCESSORIES

| | | |
|-----------------|-------------------------|-----------|
| 11-0029-9-06000 | Slide shield | Universal |
| 11-0398-9-00000 | Full extension stop | Universal |
| 11-0388-0-00000 | Pneumatic chamois liner | Custom |
| 11-0177-0-00000 | Chamois liner | Custom |
| 11-0757-x-06000 | Chamois strap pads | S, M, L |
| 11-0016-x-06000 | Sports cover | XS-XXL |
| 11-0103-x-06000 | Sports cover (short) | XS-XXL |
| 11-0075-x-00000 | Lyca undergarment | XS-XXL |
| 11-0122-x-00000 | Neoprene undergarment | XS-XXL |



NEVER STOP GETTING BETTER™

www.djortho.com