

Deluxe Elastic Ankle

Ideal Application

The Deluxe Elastic Ankle Support is a neoprene alternative, constructed of heavy-duty cotton elastic. It's ideal for warm weather climates or patients allergic to neoprene, for mild soft tissue injuries and arthritis.

Features and Benefits

- Heavy duty cotton elastic
- Breathable
- Posterior malleolar "J" buttresses for added compression

Ordering Information

PART NUMBER	
11-0346-x-xxxxx	
SIZES AVAILABLE: S, M, L	MEASUREMENTS
x=2=Small	8"-9" (20.25cm-23cm)
x=3=Medium	9"-9 3/4" (23cm-24.75cm)
x=4=Large	9 3/4"-10 1/2" (24.75cm-26.5cm)
Measure circumference directly above malleolus.	
Colors Available: white, beige	



NEVER STOP GETTING BETTER™

www.djortho.com