

## Deluxe Knee Immobilizer **Ideal Application**

The Deluxe Knee Immobilizer provides superb lightweight immobilization following injury or surgery. It is available in a wide array of sizes for excellent comfort and fit.

### Features and Benefits

- 6-strap fastening system with extra wide hook and loop
- 3 anatomically contoured aluminum posterior stays
- Anatomically contoured aluminum stay - medial and lateral



### Ordering Information

| <i>PART NUMBER</i>                     | <i>LENGTH</i>          |
|--|------------------------|
| 11-0226-x                              | 12"                    |
| 11-0227-x                              | 16"                    |
| 11-0228-x                              | 20"                    |
| 11-0229-x                              | 23"                    |
| 11-0230-x                              | 26"                    |
| Measurement-Taken 4" above mid-patella |                        |
| <i>SIZES</i>                           |                        |
| x=2=Small                              | 16" max. circumference |
| x=3=Medium                             | 19" max. circumference |
| x=4=Large                              | 22" max. circumference |
| x=5=XLarge                             | 26" max. circumference |



**NEVER STOP GETTING BETTER™**

[www.djortho.com](http://www.djortho.com)