

## Drytex Knee Support

### Ideal Application

The Drytex Knee Support sleeve is a neoprene alternative, constructed of a unique nylon core and polyester Lycra® fabric that allows for improved airflow, maintains warmth, provides compression and support yet allows for breathability. Ideal for patients who have allergies or irritation to neoprene or live in warm weather climates.

### Features and Benefits

- Available with patella or popliteal cut out
- Durable, lightweight fabric
- Alternative to neoprene

### Ordering Information

*PART NUMBER*

11-0548-x-06000	Basic Knee Sleeve
11-0658-x-06000	Basic Knee Sleeve w/open patella
11-0664-x-06060	Basic Knee Sleeve w/popliteal cut-out & open patella

*SIZES AVAILABLE: XS-XXXL Measurements taken 6" above mid-patella*

x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
x=7=XXXL= 29-1/2"-32-1/2" (75.82-82.5cm)	

Colors Available: black

Suggested HCFA Code: L1825



NEVER STOP GETTING BETTER™

[www.djortho.com](http://www.djortho.com)