

Drytex® Economy Hinged Knee

Ideal Application

For mild medial/lateral support of the knee. The Drytex Economy Hinged Knee is a neoprene alternative constructed of a unique nylon core and polyester Lycra® fabric that allows for improved airflow. This fabric provides compression and support, yet allows for breathability. Ideal for patients who have allergies or irritations to neoprene or live in warm weather climates.

Features and Benefits

- Removable polycentric aluminum uprights
- Alternative to neoprene
- Durable lightweight
- Available in wraparound style

Ordering Information

PART NUMBER	
11-0670-x-06060	Drytex Economy Hinged Knee
11-0671-x-06060	Drytex Economy Hinged Knee w/popliteal cut out
11-0672-x-06000	Drytex Economy Hinged Knee Wraparound
11-0673-x-06000	Drytex Economy Hinged Knee Wraparound w/popliteal cut out
SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella	
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
x=7=XXXL= 29-1/2"-32-1/2" (75.82-82.5cm) AVAILABLE IN WRAPAROUND ONLY	
Colors Available: black	



NEVER STOP GETTING BETTER™

www.djortho.com