

## Drytex® Hinged Knee

### Ideal Application

The Drytex Hinged Knee, designed for mild to moderate medial lateral control, is a neoprene alternative constructed of a unique nylon core and polyester Lycra® fabric that allows for improved airflow. This fabric provides compression and support, yet allows for breathability. Ideal for patients who have allergies or irritations to neoprene or live in warm weather climates.

### Features and Benefits

- Polycentric hinges
- Durable, lightweight fabric
- Extension stops: 10°, 20°, 30° and 40° (10° installed)
- Flexion stops: 45°, 60°, 75° and 90°
- Available in wraparound style

### Ordering Information

<i>PART NUMBER</i>	
11-0555-x-06000	Drytex Hinged Knee
11-0556-x-06000	Drytex Hinged Knee w/popliteal cut out
11-0669-x-06060	Drytex Hinged Knee wraparound
11-0758-x-06060	Drytex Hinged Knee wraparound w/popliteal cut out
<i>SIZES AVAILABLE: XS-XXXL Measurements taken 6" above mid-patella</i>	
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
x=7=XXXL= 29-1/2"-32-1/2" (75.82-82.5cm)	
Colors Available: black	



NEVER STOP GETTING BETTER™

[www.djortho.com](http://www.djortho.com)