

Drytex®  
Wraparound  
Playmaker®

### Ideal Application

The Drytex Wraparound Playmaker is suitable for mild to moderate ACL, PCL, MCL and LCL instabilities, and features 4-point dynamic design with ACL/PCL and CI strapping options. It is a neoprene alternative constructed of a unique nylon core and polyester Lycra® fabric that allows for improved airflow. This fabric is ideal for patients who have neoprene-based irritations or allergies, or live in warm weather climates.

### Features and Benefits

- Polycentric hinges with contoured condyle shells
- 4-point dynamic leverage system
- Extension stops: 10°, 20°, 30° and 40° (10° installed)
- Flexion stops: 45°, 60°, 75° and 90°
- Alternative to neoprene
- Lightweight, breathable fabric
- Extension/Flexion stop kit included



### Ordering Information

<i>PART NUMBER</i>	
11-0666-x-06060	
11-0759-x-06060	w/popliteal cut out
<i>SIZES AVAILABLE: XS-XXXL Measurements taken 6" above mid-patella</i>	
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
x=7=XXXL=29-1/2"-32-1/2" (75-82.5cm)	
Colors Available: black	



NEVER STOP GETTING BETTER™

[www.djortho.com](http://www.djortho.com)