

Dura*Kold® Surgical Knee Sleeve

Ideal Application

Provides extended cold therapy to the knee area for post-surgical applications. The DuraKold products can effectively penetrate surgical dressings to decrease pain, edema, and secondary hypoxic tissue injury.

Features and Benefits

- Provides 2 hours of cold therapy
- Removable ice mat inserts that can be exchanged without removing the wrap to provide continuous cold therapy
- Keeps patients comfortable and dry while allowing mobility



Ordering Information

<i>PART NUMBER</i>	
11-0904-0-02000	Surgical Knee Sleeve w/2 ice inserts
11-0905-0-02000	Surgical Knee Sleeve w/4 ice inserts
<i>ACCESSORIES</i>	
11-0888-0-02000	DuraKold Pack 9x11
11-0906-0-00000	DuraKold removable ice inserts 24/pack



NEVER STOP GETTING BETTER™

www.djortho.com