

Economy Hinged Knee with Buttress

Ideal Application

The polycentric hinge with aluminum uprights in this sleeve provide light support to medial/lateral ligaments, in addition to the support offered by neoprene and strapping.

Features and Benefits

- Polycentric hinge with 10-3/4" aluminum uprights
- 1/4" neoprene only
- Removable uprights for washing
- Available with inferior or superior buttress



Ordering Information

PART NUMBER

11-0378-x-xxxx w/popliteal cut-out & inferior buttress

SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella

x=1=XS=13"-15 1/2" (33-39cm) x=2=S=15 1/2"-18 1/2" (39-47cm)

x=3=M=18 1/2"-21" (47-53cm) x=4=L=21"-23 1/2" (53-60cm)

x=5=XL=23 1/2"-26 1/2" (60-67cm) x=6=XXL=26 1/2"-29 1/2" (67-75cm)

Colors Available: black



NEVER STOP GETTING BETTER™

www.djortho.com