

## Elastic Wrist Splint

### Ideal Application

The Elastic Wrist Splint is ideal for Carpal Tunnel Syndrome, wrist sprains and cock-up wrist immobilization.

### Features and Benefits

- Cotton-elastic construction
- Hook-and-loop closures
- Malleable aluminum cock-up splint

### Ordering Information

PART NUMBER	
11-0222-x	Right
11-0223-x	Left
SIZES AVAILABLE: S-XL MEASURE WRIST CIRCUMFERENCE	
x=2=Small=5 1/4-6" (13.25-15.25 cm)	x=4=Large 7 1/2-8 1/2" (19-21.5 cm)
x=3=Medium=6-7 1/2" (15.25-19 cm)	x=5=XLarge 8 1/2"+ (21.5 cm +)
Color Available: beige Only	

