

DonJoy® Elbow Guard

Ideal Application

A lightweight soft good with rigid elbow protection. Designed by medical professionals with your needs in mind. Provides elbow stability for mild to moderate medial/lateral and hyper-extension injuries.

Features and Benefits

- Wrap-around bicep/tricep portion with elastic strap for ease of application
- Bi-lateral hinges for flexion/extension and medial/lateral control
- "X" strap for great hyperextension control
- Bicep and tricep binding for clean finish and durability
- Elbow pad for protection in contact sports

Ordering Information

PART NUMBER	
11-1003-X-06000	DonJoy Elbow Guard
SIZES AVAILABLE: S-XXL	

*THIS BRACE CAN BE USED FOR BOTH RIGHT AND LEFT ARMS.
MEASUREMENT IS CIRCUMFERENCE OF ARM TAKEN AT MID-BICEPS.*

x=2=Small	10" to 11.5" (25.4cm to 29.21cm)
X=3=Medium	11.6" to 13" (29.464cm to 33.02cm)
x=4=Large	13.1" to 14.5" (33.274cm to 36.83cm)
x=5=XLarge	14.6" to 16" (37.084 cm to 40.64cm)
x=6=XXLarge	16.1" to 17.5" (40.894cm to 44.45cm)
x=7=XXXLarge	17.6" to 19" (44.704cm to 48.26cm)

Color Available: black



NEVER STOP GETTING BETTER™

www.djortho.com