

Telescoping ELS®

Ideal Application

The ELS features a slide lock mechanism providing the ability to lock at full extension or have free range of motion. The brace is very easy to use with a simple yet effective design. The brace is ideal for post-op use when full extension weight-bearing protocols are used. The Cool ELS gives improved patient comfort with abbreviated foam.

Features and Benefits

- Lock at full extension
- Easily size and custom fit the brace to the patient
- Brace can extend to ankle to help prevent migration
- 50 degree flexion stop kit available for meniscal repair protocols

Ordering Information

<i>PART NUMBER</i>	<i>SIZE</i>
11-0965-9-13066 Enhanced ELS, Standard	
Color Available: grey	
<i>ACCESSORIES</i>	
11-0080-9-00000	50 flexion stop, ELS
11-0024-x	Ankle stirrup, right
11-0026-x	Ankle stirrup, left



NEVER STOP GETTING BETTER™

www.djortho.com