

Female Fource - Custom

Ideal Application

Females are 3-8 times more likely to sustain an ACL tear than males. In addition, over 70 percent of all ACL injuries occur during non-contact activities, such as cutting, jumping and twisting maneuvers. The DonJoy® Female Fource™ addresses biomechanical factors that have been implicated in the high incidence of anterior cruciate ligament (ACL) tears in females. Females tend to perform cutting and landing maneuvers with their knees and hips in an extended and less protected position. The Female Fource incorporates dj Orthopedics' FourcePoint™ hinge technology, which has been clinically demonstrated to increase knee flexion angles, limiting the time the knee spends in at-risk positions during activity; therefore, providing females the protection they need.

The Female Fource is the first functional knee brace to come with a custom graphics sticker sheet to allow patients to personalize their braces to their own liking. In addition, with nearly 212,000 females being diagnosed with breast cancer and approximately 40,000 losing their life from the disease in 2005, DonJoy will be contributing to breast cancer initiatives.

The Custom Defiance III and Defiance Female Fource is DonJoy's premier knee brace for females. Combining the ultimate in technology, materials and design, this lightweight brace provides optimal support for moderate to severe female ACL instabilities, ACL reconstructions, collateral ligament instabilities, hyperextension and for prophylactic use. The low-profile design allows for a comfortable fit under clothing. As the strongest of DonJoy's functional braces, the custom, carbon composite Defiance III and Defiance Female Fource is recommended for high collision sports, while its featherweight materials and custom fit also make it the brace of choice for activities of daily living.



Features and Benefits

- Patented FourcePoint Hinge Technology: The only bracing technology clinically proven to protect the ACL*
 - Reduces the time spent near full extension, or in the "at risk" position
 - Enhances the Four-Points-of-Leverage design by increasing the posterior load applied by the brace on the tibia
 - Improves a patient's natural gait by applying progressive resistance before reaching a firm end point
- Four-Points-of-Leverage design for superior ligament protection and comfort
- Medial Swooping Thigh Frame (Defiance III) for increased medial thigh clearance
- Swiveling strap tabs (Defiance III) to match anatomy through movement
- Anti-Migration calf strap and Supra condylar Pad for brace suspension
- Custom Graphics Sticker Kit to personalize the brace to your own liking
- Lightweight (18oz), low profile design
- Short calf versions available
- Short thigh available (Defiance)
- Assortment of color, graphic and gradation options available
- Lifetime limited warranty on frame and hinge
- 24 Hour turnaround time
- Breast Cancer Donation to support breast cancer initiatives

Ordering Information

ACCESSORIES:

11-0075-x-00000	Lycra Undergarment
11-0122-x	Neoprene Undergarment
11-0016-x-06000	Sports Cover, Standard
11-0103-x-06000	Sports Cover, Short
11-2005-x	Anti-Migration Calf Band
11-0543-x-06000	Neoprene Suspension Strap Kit
11-1367-9-06000	Calf Pinch Guard
25-1923-0-06000	Air Condyle Kit (Qty.2)
25-1480-0-00000	Brace Bag

SIZING INFORMATION	THIGH	CALF
EACH CUSTOM MADE.	MINIMUM 12 1/2" (31.8CM)	10 1/2" (26CM)
EACH CUSTOM MADE.	MAXIMUM 38 1/8" (96.83CM)	26 3/4" (67.94 CM)

Custom braces are available in an assortment of colors and designs



NEVER STOP GETTING BETTER™

www.djortho.com