

Four Point

Ideal Application

Fource Point Technology was developed to further protect the ACL and enhance the Four Points of Leverage brace design. The patented Fource Point hinges are designed to react during leg extension and offer an increasing resistance in the last 35 degrees of extension, in order to reduce extension shock and increase posterior tibial loads. The hinge has three adjustment settings to vary the amount of resistance to meet the specification of the physician's protocols and the patient's comfort.

The Fource Point is ideal for the patient with ACL reconstruction, ACL deficiencies or prophylactic use.



Features and Benefits

- Adjustable resistance offered in 25°-0°, 35°-10°, 45°-20° (shipped with 35°-10°)
- Gradual "soft" stop to extension
- Enhanced Four Points of Leverage design
- Standard anti migration Supra Condyle pad
- 7 sizes and three lengths
- Moldable aluminum cuffs for custom fit
- Removable thermal-formed liners and strap pads
- Flexion/extension stops
- One-year warranty on frame and hinge / six months on softgoods

Ordering Information

PART NUMBER

11-1300-x-06060 (1-7)	Fource Point Right ACL Standard
11-1301-x-06060 (1-7)	Fource Point Left ACL Standard
11-1302-x-06060 (1-7)	Fource Point Right ACL Short
11-1303-x-06060 (1-7)	Fource Point Left ACL Short

Colors Available:black

SIZING INFORMATION: XS-XXXL

x=1=XS	13"-15 1/2" (33-39cm)
x=2=S	15 1/2"-18 1/2" (39-47cm)
x=3=M	18 1/2"-21" (47-53cm)
x=4=L	21"-23 1/2" (53-60cm)
x=5=XL	23 1/2"-26 1/2" (60-67cm)
x=6=XXL	26 1/2"-29 1/2" (67-75cm)
x=7=XXXL	29 1/2"-32" (75-81.25 cm)

MEASURING INFORMATION - Measurements taken 6" above mid-patella

ACCESSORIES

11-0194-9-06000	Patella cup
11-0398-9-00000	Full extension stop
11-0103-x-06000	Sports cover (short)
11-0016-x-06000	Sports Cover (standard)
11-0075-x-00000	Lycra undergarment



NEVER STOP GETTING BETTER™

www.djortho.com