

Hinged Air DonJoy® Ideal Application

The Drytex Hinged Air DonJoy is designed to treat and ease pain associated with patellofemoral malalignment, chondromalacia and tendonitis plus medial/lateral instabilities and post-operative meniscectomies. The Hinged Air DonJoy is constructed with Drytex, a neoprene alternative material that allows for improved air flow.

Features and Benefits

- Inflatable pneumatic buttresses for medial, lateral or full patellar support
- Alternative to tape and neoprene
- Drytex material for added patient comfort
- Half wrap on posterior thigh for ease of application
- Popliteal cutout
- Polycentric hinge with flexion stops of 0°, 10°, 20°, 30° and 40° (10° installed)
- Extension stops 45°, 60°, 75°, 90°

Ordering Information

<i>PART NUMBER</i>	
11-0855-x-06060	
<i>SIZES AVAILABLE: XS-XXL</i>	
<i>Measurements taken 6"</i>	
<i>above mid-patella</i>	
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
Colors Available: black	



NEVER STOP GETTING BETTER™

www.djortho.com