

Hinged "H" Buttress

Ideal Application

Designed for medial/lateral instabilities with patella tendonitis, subluxation or chondromalacia. The unique sewn-in "H" buttress provides compressive inferior and medial/lateral patellar support. Anterior "bib" completely covers open patella and hinges.

Features and Benefits

- Available with open pop only
- 1/8" neoprene
- Sewn-on anterior "bib" covers polycentric hinges
- Flexion stops of 10°, 20°, 30°, 40°, (10° installed)
- Extension stops of 45°, 60°, 75°, 90°

Ordering Information

<i>PART NUMBER</i>	
11-0828-x-06060	
<i>SIZES AVAILABLE: XS-XXXL Measurements taken 6" above mid-patella</i>	
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
x=7=XXXL=29 1/2"-32" (75-81.25cm)	
Colors Available: black	



NEVER STOP GETTING BETTER™

www.djortho.com