

Telescoping IROM®

Ideal Application

The IROM is the industry standard providing precise range-of-motion control in 10 degree increments. The IROM is ideal for ACL, PCL, MCL, and LCL surgeries. The telescoping bars allow the brace to fit a variety of patients.

Features and Benefits

- Precise range of motion control
- Easily size and custom fit the brace to the patient
- Buckles allow for fast and easy strapping
- Brace can extend to ankle to help prevent migration



Ordering Information

<i>PART NUMBER</i>	<i>SIZE</i>
11-0968-9-13066	Enhanced IROM, Standard
11-0971-9-13066	Enhanced IROM, XL XL fits 26-32" circumference/girth

Color available: grey

<i>ACCESSORIES</i>	
11-0024-x	Ankle stirrup, right
11-0026-x	Ankle stirrup, left



NEVER STOP GETTING BETTER™

www.djortho.com