

IROM® Bilateral with and w/out shells

Ideal Application

For added stiffness and tissue containment, the IROM with shells offers full-circumference shells. The knee brace offers effective post-operative immobilization and protected range of motion associated with ligament and fracture related surgeries. The IROM Bilateral is a short 18-inch version and is available with and without shells.

Features and Benefits

- 10 degree incremental range of motion control
- Quick-release buckles allow for fast and easy strapping
- Shells provide added stability



Ordering Information

PART NUMBER

11-2127-9-13066	Enhanced IROM Bilateral
11-0972-9-13066	Enhanced IROM Bilateral w/ shells



NEVER STOP GETTING BETTER™

www.djortho.com