

Knee Support Ideal Application

Constructed of high-quality neoprene or Drytex for compression and support.

Features and Benefits

- 1/4" or 1/8" neoprene or Drytex material for improved airflow and compressive support
- Available with patella or popliteal cut-out



Ordering Information

PART NUMBER	PART NUMBER	
1/4" NEOPRENE	1/8" NEOPRENE	
11-0008-x-xxxxx	11-0002-x-xxxxx	double-sided nylon
11-0108-x	11-0102-x	single-sided nylon (blue only)
11-0358-x-xxxxx	11-0058-x-xxxxx	double-sided nylon w/patella cut-out
11-0370-x-xxxxx	11-0368-x-xxxxx	double sided nylon w/popliteal cut-out
	11-0371-x-xxxxx	w/patella and popliteal cut-outs
11-0548-x-06000	Drytex® Basic Knee Sleeve	
11-0658-x-06000	Drytex® Basic Knee Sleeve w/open patella	
11-0664-x-06000	Drytex® Basic Knee Sleeve w/popliteal cut-out & open patella	
<i>SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella</i>		
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)	
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)	
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)	
Colors Available: blue, black		

Suggested HCFA Code: L1825



NEVER STOP GETTING BETTER™

www.djortho.com