

Legend™

Ideal Application

The Legend™ by DonJoy has made a name for itself as a strong, low-profile, custom-fit brace. Whether your passion is contact sports, salt- or fresh-water sports activities, or simply to get the most out of the activities of daily living, the Legend™ can take whatever you give. Excellent for those with moderate to severe ACL, PCL, MCL, and LCL instabilities, the Legend™ by DonJoy will support your knees and support your goal to Never Stop Getting Better.

Features and Benefits

- ACL and PCL versions feature patented 4-point dynamic leverage system
- Low-profile design
- CI version available
- 7 sizes (ranging from a 13" to 32" thigh)
- Removable doeskin liner (extra set included)
- Pneumatic condyle pads for added comfort and fit
- Extension stops: 0°, 10°, 20°, 30° and 40° (10° installed)
- Flexion stops: 45°, 60°, 75° and 90°
- Chamois liners and strap pads available
- One-year warranty on frame and hinge

Ordering Information

| <i>RIGHT PART NUMBER</i> | <i>VERSION</i> | <i>LEFT PART NUMBER</i> | <i>VERSION</i> |
|------------------------------|----------------|-----------------------------|----------------|
| 11-0422-x-06060 | ACL | 11-0423-x-06060 | ACL |
| 11-0442-x-06060 | PCL | 11-0443-x-06060 | PCL |
| 11-0452-x-06060 | CI | 11-0453-x-06060 | CI |

| <i>SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella</i> | |
|--|-----------------------------------|
| x=1=XS=13"-15 1/2" (33-39cm) | x=2=S=15 1/2"-18 1/2" (39-47cm) |
| x=3=M=18 1/2"-21" (47-53cm) | x=4=L=21"-23 1/2" (53-60cm) |
| x=5=XL=23 1/2"-26 1/2" (60-67cm) | x=6=XXL=26 1/2"-29 1/2" (67-75cm) |
| x=7=XXXL=29 1/2"-32" (75-81.25 cm) | |

Colors available: black

ACCESSORIES

| | | | |
|---------------------------------------|--------|-------------------------------------|------|
| 11-0194-9-06000 Patella cup | Univ | 11-0398-9-00000 Full extension stop | Univ |
| 11-0016-x-06000 Sports cover | XS-XXL | | |
| 11-0075-x-00000 Lycra undergarment | XS-XXL | | |
| 11-0122-x-00000 Neoprene undergarment | XS-XXL | | |



NEVER STOP GETTING BETTER™

www.djortho.com