

## Lerman Shoulder Ideal Application

For immobilization following shoulder fusion, dislocation, reconstruction or soft tissue repairs. This brace allows the shoulder to be positioned in external, neutral or internal rotation. Following brachial plexus injury this orthosis can be used to prevent further stretch during the healing process and can position the hand for functional activities.

### Features and Benefits

- Fully adjustable in all planes - allows for desired specific abduction up to 45°
- Comfortable
- Easy to apply
- One size fits most (28"-40")
- Universal - fits left or right
- Patented open cell non-slip material for patient comfort and consistent position



### Ordering Information

<i>PART NUMBER</i>		
11-1361-9	Lerman Shoulder Orthosis	Universal Size



NEVER STOP GETTING BETTER™

[www.djortho.com](http://www.djortho.com)