

Lumbosacral Support

Ideal Application

Ideal for lumbar support, mild to moderate low back pain, and low back sprains or strains as well as abdominal compression.

Features and Benefits

- Semi-rigid low back support
- 4-heat formable thermo-plastic stays for individual custom fit
- Warm and compressive neoprene
- Additional side pulls for added support and compression



Ordering Information

PART NUMBER

11-0015-x-xxxxx

SIZES AVAILABLE: S-XXL

x=2=Small=24-30" (61-75 cm)

x=3=Medium=30-36" (75-90 cm)

x=4=Large=36-42" (90-105 cm)

x=5-XLarge=42-48" (105-120 cm)

x=6=XXLarge=48-54" (120-135 cm)

Colors Available: blue, black

Suggested HCFA Code: L0515



NEVER STOP GETTING BETTER™

www.djortho.com