

MaxTrax Walker Ideal Application

Designed with a wider foot bed for comfort and stability. A low profile, rocker bottom helps promote natural gait and reduce plantar pressures. Cushioned inner and outer sole provides shock absorption to help improve patient comfort during ambulation. Malleable uprights breakdown on standard height.

Features and Benefits

- Wider footbed for increased stability
- Swiveling D-rings allow for criss cross or parallel strapping
- Cushioned inner/outer sole are designed to absorb shock upon heel strike
- Breakdown uprights for short tibia/fibula
- One-piece design keeps seams away from post surgical areas
- Anterior/posterior and medial/ lateral shells for compression and support

Ordering Information

PART NUMBER	PRODUCT DESCRIPTION	SIZE
11-1382-x-00000	MaxTrax Walker	1 - 5

Colors available: black

ACCESSORIES

PART NUMBER	PRODUCT DESCRIPTION	SIZE
11-1372-x-00000	Liner Pad Kit	1 - 5
11-1373-x-00000	Replacement Liner Kit, Tall Walkers	1 - 5
11-1377-x-00000	Replacement Liner Kit, Short Walkers	1 - 5
25-3772-x-06000	Replacement Air Liner/ Tall Walkers	1 - 5
25-3773-x-06000	Replacement Air Liner, Short Walkers	1 - 5
11-1375-x-00000	Anterior/Posterior & Medial/Lateral Shells	2 = XS - S, 5 = M - XL
11-1374-x-00000	Strap Replacement Kit	1 - 5

SIZING INFORMATION

PRODUCT SIZE	MEN SHOE SIZE	WOMEN SHOE SIZE	PEDIATRIC
X=1=X Small	NA	Up to 4	Up to 2
X=2=Small	up to 5	4.5 - 6	
X=3=Medium	5.5 - 10	6.5 - 11	
X=4=Large	10.5 - 13.5	11.5 - 14.5	
X=5=X Large	14 - 17	NA	

Measuring Information: Measurements 6" above mid patella



NEVER STOP GETTING BETTER™

www.djortho.com