

## MaxTrax Air Ankle Walker

### Ideal Application

The MaxTrax Air Ankle walkers are designed to increase stability and help decrease pain and edema following trauma or post-operative procedures. Adjustable air chambers to accommodate changes in edema during rehab and maximize patient compliance.

### Features and Benefits

- Adjustable air chambers to accommodate changes in edema during rehab
- Maximize patient compliance
- Integrated pump designed to make the liner simple to inflate and adjust
- Medial/lateral dual air cell system inflates sides equally

### Ordering Information

PART NUMBER	PRODUCT DESCRIPTION	SIZE
11-1380-x	MaxTrax Air Ankle Walker	1 - 5

**Colors available: black**

**ACCESSORIES**

PART NUMBER	PRODUCT DESCRIPTION	SIZE
11-1372-x-00000	Liner Pad Kit	1 - 5
11-1373-x-00000	Replacement Liner Kit, Tall Walkers	1 - 5
11-1377-x-00000	Replacement Liner Kit, Short Walkers	1 - 5
25-3772-x-06000	Replacement Air Liner/ Tall Walkers	1 - 5
25-3773-x-06000	Replacement Air Liner, Short Walkers	1 - 5
11-1375-x-00000	Anterior/Posterior & Medial/Lateral Shells	2 = XS - S, 5 = M - XL
11-1374-x-00000	Strap Replacement Kit	1 - 5

**SIZING INFORMATION**

PRODUCT SIZE	MEN SHOE SIZE	WOMEN SHOE SIZE	PEDIATRIC
X=1=X Small	NA	Up to 4	Up to 2
X=2=Small	up to 5	4.5 - 6	
X=3=Medium	5.5 - 10	6.5 - 11	
X=4=Large	10.5 - 13.5	11.5 - 14.5	
X=5=X Large	14 - 17	NA	

**Measuring Information: Measurements 6" above mid patella**



NEVER STOP GETTING BETTER™

[www.djortho.com](http://www.djortho.com)