

MaxTrax ROM Walker

Ideal Application

Provides protected range-of-motion following trauma or post-operative procedures. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights will lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar- and dorsi-flexion.

Features and Benefits

- Pediatric sizing through extra large
- Wider footbed for increased stability
- Swiveling D-rings allow for criss cross or parallel strapping
- Cushioned inner/outer sole are designed to absorb shock upon heel strike
- Breakdown uprights for short tibia/fibula
- One-piece design keeps seams away from post surgical areas
- Anterior/posterior and medial/ lateral shells for compression and support

Ordering Information

PART NUMBER	PRODUCT DESCRIPTION	SIZE
11-1383-x-00000	MaxTrax ROM Walker	1-5

Colors available: black

ACCESSORIES

PART NUMBER	PRODUCT DESCRIPTION	SIZE
11-1372-x-00000	Liner Pad Kit	1 - 5
11-1373-x-00000	Replacement Liner Kit, Tall Walkers	1 - 5
11-1377-x-00000	Replacement Liner Kit, Short Walkers	1 - 5
25-3772-x-06000	Replacement Air Liner/ Tall Walkers	1 - 5
25-3773-x-06000	Replacement Air Liner, Short Walkers	1 - 5
11-1375-x-00000	Anterior/Posterior & Medial/Lateral Shells	2 = XS - S, 5 = M - XL
11-1374-x-00000	Strap Replacement Kit	1 - 5

SIZING INFORMATION

PRODUCT SIZE	MEN SHOE SIZE	WOMEN SHOE SIZE	PEDIATRIC
X=1=X Small	NA	Up to 4	Up to 2
X=2=Small	up to 5	4.5 - 6	
X=3=Medium	5.5 - 10	6.5 - 11	
X=4=Large	10.5 - 13.5	11.5 - 14.5	
X=5=X Large	14 - 17	NA	

Measuring Information: Measurements 6" above mid patella

