

MaxTrax ROM Air Walker

Ideal Application

Provides protected range-of-motion following trauma or post-operative procedures. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights will lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar- and dorsi-flexion.

Features and Benefits

- Adjustable air chambers to accommodate changes in edema during rehab
- Maximize patient compliance
- Integrated pump designed to make the liner simple to inflate and adjust
- Medial/lateral dual air cell system inflates sides equally

Ordering Information

| PART NUMBER | PRODUCT DESCRIPTION | SIZE |
|-------------|------------------------|------|
| 11-1381-x | MaxTrax ROM Air Walker | 1-5 |

Colors available: black

ACCESSORIES

| PART NUMBER | PRODUCT DESCRIPTION | SIZE |
|-----------------|--|------------------------|
| 11-1372-x-00000 | Liner Pad Kit | 1 - 5 |
| 11-1373-x-00000 | Replacement Liner Kit, Tall Walkers | 1 - 5 |
| 11-1377-x-00000 | Replacement Liner Kit, Short Walkers | 1 - 5 |
| 25-3772-x-06000 | Replacement Air Liner/ Tall Walkers | 1 - 5 |
| 25-3773-x-06000 | Replacement Air Liner, Short Walkers | 1 - 5 |
| 11-1375-x-00000 | Anterior/Posterior & Medial/Lateral Shells | 2 = XS - S, 5 = M - XL |
| 11-1374-x-00000 | Strap Replacement Kit | 1 - 5 |

SIZING INFORMATION

| PRODUCT SIZE | MEN SHOE SIZE | WOMEN SHOE SIZE | PEDIATRIC |
|--------------|---------------|-----------------|-----------|
| X=1=X Small | NA | Up to 4 | Up to 2 |
| X=2=Small | up to 5 | 4.5 - 6 | |
| X=3=Medium | 5.5 - 10 | 6.5 - 11 | |
| X=4=Large | 10.5 - 13.5 | 11.5 - 14.5 | |
| X=5=X Large | 14 - 17 | NA | |

Measuring Information: Measurements 6" above mid patella

