

MC Walker®

Ideal Application

The Motion Control (MC) Walker is recommended for use following Achilles tendon repairs, stable fractures of the lower leg and ankle, acute ankle sprains, and soft tissue injury. It provides foot and ankle immobilization in 0°, 10° or 20° of plantarflexion as well as protected range of motion in 7.5° increments from 45° plantarflexion to 30° dorsiflexion. Extension pads for the liner are available to accommodate larger patients.

Features and Benefits

- Rocker bottom sole
- Washable liner
- 7.5° range-of-motion increments
- Low profile
- Promotes early range-of-motion



Ordering Information

<i>PART NUMBER</i>		
11-0459-x-06136		
<i>SIZES AVAILABLE: S,M,L</i>		
	<i>MEN</i>	<i>WOMEN</i>
x=2=Small=Shoe size:	up to 6	up to 7
x=3=Medium=Shoe size:	6-10	7-11
x=4=Large=Shoe size:	10+	11+
<i>ACCESSORIES</i>		
11-0825-9-13000	Extension pad	6-pack
<i>Suggested HCFA Code: L2114 + L2200 or L1930 + L2200</i>		



NEVER STOP GETTING BETTER™

www.djortho.com