

Playmaker®

Ideal Application

Suitable for mild to moderate ACL and/or PCL, MCL and LCL instabilities, this neoprene brace with extended hinge bars is an ideal intermediate-level product that can be used in sporting activities.

Features and Benefits

- ACL and PCL versions feature patented 4-point dynamic leverage system
- Low profile
- Polycentric hinge with contoured condyle shells
- Versatile strapping
- Extension stops: 10°, 20°, 30° and 40° (10° installed)
- Flexion stops: 45°, 60°, 75° and 90°
- Additional full-circumference thigh strap for CI conversion



Ordering Information

PART NUMBER

11-0860-x-xxxxx

11-0861-x-xxxxx w/popliteal cut out

SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella

x=1=XS=13"-15 1/2" (33-39cm)

x=2=S=15 1/2"-18 1/2" (39-47cm)

x=3=M=18 1/2"-21" (47-53cm)

x=4=L=21"-23 1/2" (53-60cm)

x=5=XL=23 1/2"-26 1/2" (60-67cm)

x=6=XXL=26 1/2"-29 1/2" (67-75cm)

Colors Available: black

ACCESSORIES

12-0422-0-00000

Full extension stop kit



NEVER STOP GETTING BETTER™

www.djortho.com