

Playmaker® with Patella Donut

Ideal Application

The Playmaker is designed for mild ACL, PCL or Combined Instabilities (CI) with patella tendonitis, chondromalacia, or subluxation. The supportive donut may be trimmed to allow for inferior, superior, medial or lateral buttress.

Features and Benefits

- Patented 4-point configuration
- Available with or without popliteal cutout
- Trimable donut support
- Flexion stops of 0°, 10°, 20°, 30°, 40° (10° installed)
- Extension stops of 45°, 60°, 75°, 90°

Ordering Information

PART NUMBER	
11-0862-x-xxx66	w/donut
11-0863-x-xxx66	w/donut and popliteal cutout
SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella	
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
Colors Available: black	



NEVER STOP GETTING BETTER™

www.djortho.com