

Playmaker®
Wraparound/
Wraparound
IROM
Playmaker®

Ideal Application

The wraparound design makes this brace ideal for those patients who have trouble applying a sleeve-style brace. It provides dependable support for mild to moderate ACL and/or PCL, MCL and LCL instabilities. The IROM Playmaker is ideal for post-op care following meniscal repairs as well as ACL/PCL reconstruction.

Features and Benefits

- ACL and PCL versions feature patented 4-point dynamic design
- Wraparound style neoprene sleeve
- Available with polycentric or IROM hinge
- Polycentric hinge provides extension stops of 10°, 20°, 30° and 40° (10° installed) and flexion stops of 45°, 60°, 75° and 90°
- Additional full-circumference thigh strap for CI conversion

Ordering Information

<i>PART NUMBER</i>	
11-0864-x-xxxxx	
11-0865-x-xxxxx	w/popliteal cut out
<i>WRAPAROUND IROM PLAYMAKER®</i>	
11-0866-x-xxxxx	
11-0867-x-xxxxx	w/popliteal cut out
<i>SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella</i>	
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
Colors Available: black	



NEVER STOP GETTING BETTER™

www.djortho.com