

Post-Op Ankle Stirrup

Ideal Application

This accessory is designed to fit with DonJoy's long leg post-op braces, attaching easily inside the hinge bars to extend brace leverage for added rotational control. Contoured heel cup is tapered for comfort. Free pivoting grommets permit unrestricted dorsi/plantar flexion. A dorsum strap keeps the heel securely and comfortably seated.

Features and Benefits

- Compatible with all DonJoy long leg braces
- Contoured heel cup
- Unrestricted dorsi and plantarflexion

Ordering Information

PART NUMBER

11-0024-x Right

11-0026-x Left

SIZES AVAILABLE: S, M, L

x=2=Small=Men's shoe size up to 8

x=3=Medium=Men's shoe size 8 -10

x=4=Large=Men's shoe size 10+

