

# ROCKETSOC®

## ANKLE BRACE

### Patients

- Active adolescents
- Active adults

### Indications

- Acute sprains
- Chronic instabilities
- Prophylactic use

### Features and Benefits

- Patented proven strapping system locks calcaneus and forefoot
- Protects 3 most commonly injured ligaments
  - Anterior Talofibular Ligament
  - Posterior Talofibular Ligament
  - Calcaneofibular Ligament
- Right and Left versions
- Durability for all sports
- Fits easily into an athletic shoe
- Adjustable strapping allows maximum support and comfort
- Medial/Lateral superior malleolus buttress for comfort and suspension



Recommended HCPCS Code: L1906



## Strap In. Take Off.

The RocketSoc is DonJoy's most durable and supportive soft ankle brace. The unique patented strapping system locks the ankle to prevent lateral ankle sprains without restricting range of motion.

### LACE UP



- Lace up design allows contouring fit
- Low profile lightweight Nylon construction
- Fewest straps of all RocketSoc models

### Ordering Information

Part Number	Ankle	Color
11-0338-x-06000	Right	Black
11-0339-x-06000	Left	Black

### NEOPRENE



- Wraparound style
- Neoprene fabric for even compression and warmth
- Optional inferior buttress for swelling

### Ordering Information

Part Number	Ankle	Color
11-0336-x-06000	Right	Black
11-0337-x-06000	Left	Black

### DRYTEX®



- Wraparound style
- Breathable Drytex® material keeps ankle cool and dry
- Optional inferior buttress for swelling

### Ordering Information

Part Number	Ankle	Color
11-0665-x-06000	Right	Black
11-0667-x-06000	Left	Black

## SIZING

Brace Size	Men's Shoe Size	Women's Shoe Size
x = 1 = XS	4 – 5.5	5 – 7.5
x = 2 = S	6 – 7.5	8 – 9.5
x = 3 = M	8 – 9.5	10 – 11.5
x = 4 = L	10 – 11.5	12 – 13.5
x = 5 = XL	12 – 13.5	N/A
x = 6 = XXL	14 – 16.5	N/A

dj ortho® 1430 Decision Street, Vista, CA 92081-8553  
 t800.336.6569 f800.936.6569  
 www.donjoy.com © 2006 dj Orthopedics, LLC



It is within the sole discretion of the customer to determine the appropriate billing code for a product, as well as, whether the use of a product complies with medical necessity and other documentation requirements of the payer. dj Orthopedics, LLC accepts no responsibility whatsoever in this regard.  
 U.S. Patent 5,330,419 & 5,217,431