

S.C.O.I. Shoulder Brace

Ideal Application

The S.C.O.I. shoulder brace offers immobilization for rotator cuff repairs, capsular shifts, Bankhart procedures, and prosthetic shoulder replacements. The S.C.O.I. allows for 30°-150° of abduction/adduction.

Features and Benefits

- Allows 30°-150° abduction/adduction range-of-motion hinge
- Universal, fits right or left
- No auxiliary hardware
- Elbow range-of-motion hinge



Ordering Information

PART NUMBER

11-0100-9-13136

UNIVERSAL: FITS RIGHT OR LEFT



NEVER STOP GETTING BETTER™

www.djortho.com