

Spiral lastic Knee

Ideal Application

Patients allergic to neoprene can still get the same even compression with this cotton elastic sleeve. Two medial/lateral spiral stays provide additional support.

Features and Benefits

- Seamless two-way stretch cotton elastic
- Open patella
- Even compression without heat
- Two spiral stays

Ordering Information

PART NUMBER

11-0240-x

SIZES AVAILABLE: S-XL Measurements taken 6" above mid-patella

x=2=S=15 1/2"-18 1/2" (39-47cm) x=3=M=18 1/2"-21" (47-53cm)

x=4=L=21"-23 1/2" (53-60cm) x=5=XL=23 1/2"-26 1/2" (60-67cm)



NEVER STOP GETTING BETTER™

www.djortho.com