

## Tru-Pull®

### Ideal Application

A simple, easy to apply bracing solution for patellar subluxation, dislocation, malalignment and instabilities.

The Tru-Pull Wraparound brace, designed by Dr. John Fulkerson, incorporates his patented technology of independent straps, which provide a consistent pull on the patella. A choice of buttresses give the patient an option for more aggressive support. The Tru-Pull Wraparound is an ideal choice for patellar subluxation, dislocation, malalignment and instabilities.

### Features and Benefits

- Patented bifurcated pull strap design allows for constant pull on patella.
- Choice of two adjustable thickness' of buttresses (1/4" and 3/8").
- Universal left and right.
- Wrap style design for ease of application and fit.
- Non-slip neoprene.



### Ordering Information

*PART NUMBER*

11-1111-X-06000 Tru-Pull Wrap

11-1112-0-06000 Buttress Kit Universal

*SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella*

x=2=S=15 1/2"-18 1/2" (39-47cm)    x=3=M=18 1/2"-21" (47-53cm)

x=4=L=21"-23 1/2" (53-60cm)    x=5=XL=23 1/2"-26 1/2" (60-67cm)



**NEVER STOP GETTING BETTER™**

[www.djortho.com](http://www.djortho.com)