

Tru-Pull™ Advanced System

Ideal Application

For symptomatic pain relief of patellofemoral dysfunction including subluxation/dislocation, mal-tracking and patellar tendonitis.

Features and Benefits

- Elastomeric pull device, with pull straps above and below the patella, has a rubber band type effect providing a strong dynamic pull on the patella as the knee extends.
- Independent anchors prevent rotation of the brace and provide a "True" and consistent pull on the patella to avoid lateral subluxation.
- Trimable adjustable buttress allows for versatility of contact area.
- Choice of 1/4" and 3/8" buttresses to give the patient the leverage they require.
- Breathable material makes the brace light and comfortable.
- Versions with and without popliteal cutout for patient preference.
- New unique materials give it a modern nonmedical look.
- Ease of application leads to excellent patient compliance.



Ordering Information

<i>PART NUMBER</i>	<i>PRODUCT DESCRIPTION</i>
11-1289-x-06000	Tru-Pull Sleeve RT
11-1290-x-06000	Tru-Pull Sleeve LT
11-1292-x-06000	Tru-Pull Sleeve w/pop RT
11-1293-x-06000	Tru-Pull Sleeve w/pop LT
<i>COLORS AVAILABLE:</i>	
<i>black</i>	
<i>SIZING INFORMATION:</i>	
<i>x=1=XS</i>	<i>13" - 15 1/2" (33-39 cm)</i>
<i>x=2=S</i>	<i>15 1/2" - 18 1/2" (39-47 cm)</i>
<i>x=3=M</i>	<i>18 1/2" - 21" (47-53 cm)</i>
<i>x=4=L</i>	<i>21" - 23 1/2" (53-60 cm)</i>
<i>x=5=XL</i>	<i>23 1/2" - 26 1/2" (60-67 cm)</i>
<i>x=6=XXL</i>	<i>29 1/2" - 32" (75-81 cm)</i>
<i>MEASURING INFORMATION:</i>	
Measurements 6" above mid patella	
<i>ACCESSORIES:</i>	
25-3715-0-06000	1/4" Buttress
25-3716-0-06000	3/8" Buttress