

Walkabout®

Ideal Application

The Walkabout is the newest edition to the DonJoy family of walkers. This simple one piece injection molded walker is designed to be durable and comfortable yet economical. The lowest profile rocker bottom in the DonJoy line allows for normal ambulation. Ideal for acute ankle sprains, soft tissue injuries, stress fractures of the lower leg and ankle as well as stable fractures of the ankle. Extension pads for the liner are available to accommodate larger patients.

Features and Benefits

- Light weight
- Rocker bottom
- Low profile
- Economical
- Integrated straps
- Breathable fiber liner



Ordering Information

<i>PART NUMBER</i>		
11-0461-x-xxxxx		
<i>SIZES AVAILABLE: S, M, L</i>	<i>MEN'S</i>	<i>WOMEN'S</i>
x=2=S=Shoe size:	4-6	up to 7
x=3=Medium=Shoe size:	6-10	7-11
x=4=Large=Shoe size:	10-14	11+
<i>ACCESSORIES</i>		
11-0825-9-13000	Extension pad	6-pack



NEVER STOP GETTING BETTER™

www.djortho.com