

Wrist Wraps Ideal Application

The Wrist Wraps Support is ideal when limitation of range is needed in one direction only, such as extension, flexion, ulnar or radial deviation. It is also ideal for sports-related and job-related injuries allowing freedom of the hand from splinting. Can be used in football (linemen) for prophylactic use as well as jobs requiring repetitive hand/wrist movement.

Features and Benefits

- Neoprene inner lining for comfort
- Interchangeable foam inserts for variable support
- Woven nylon outer shell for durability
- Opposing straps limit rotation of brace

Ordering Information

PART NUMBER

11-0386-x-xxxxx

SIZES AVAILABLE: S, M, L UNIVERSAL: FITS RIGHT OR LEFT

x=2= Small=5-6" (12.75-15.25 cm)

x=3=Medium=6-7" (15.25-17.75 cm)

x=4=Large=7-8" (17.75 cm - 20.25 cm)

Colors Available: blue, black



NEVER STOP GETTING BETTER™

www.djortho.com